
# Tips for Qiagen® QIAgility



**Product Highlights**

- Suitable for QIAgility® (formerly CAS-1200)
- Suitable for Tecan® Freedon EVO Genesis® series and MiniPrep®

## Product Applications

PCR Setup
Assay-Setup
Serial Dilutions
Distributing Reagents
Sample or Reagent Transfer

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk

# General Data

| | |
|---|---|
| Art. No. | See variations |
| Sterile | Yes |
| Filter tip | Filter |
| Pack Size | 96 Tips (96 Tips) |
| Rack Style | Clear Hanging Rack |
| Conductibility | Yes |
| Head | Single-Channel Pipette |
| For use with (SL) | Qiagen® |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk

# More informations about Tips for Qiagen® QIAgility

STARLAB robotic tips are engineered to the same specifications as the original equipment manufacturers' tips to ensure compatibility and a consistent fit. Don't risk having your time and results compromised with defective materials and poor performance! With STARLAB, you can put your trust in disposables from a company that offers much more than just tips.

All robotic tips are certified RNase, DNase and DNA free

**Quality control**
Quality control is the essential factor which gives you the freedom to walk away from your automated workstation and get on with other tasks. You need the assurance that the tip will not let you down. STARLAB robotic tips are put through stringent and thorough testing. Every lot of tips is subject to an array of demanding tests, from the geometry of the tip itself to automation-specific requirements such as clean ejection from the liquid handling head. Every lot is required to meet or exceed the parameters necessary to give you that confidence to walk away.

Tips are tested on the original robotic equipment so there is no doubt they will perform to the same requirements as the disposables from the robot manufacturer themselves.

**No special set-up required!**
STARLAB robotic tips require no instrument hardware of software adjustment, thereby providing a simple, safe and direct high-performance substitute.

**Purchasing flexibility**
STARLAB likes to make things as simple as possible for its customers. We realise that you may need to order robotic consumables in large quantities, but may not have the space to store them. STARLAB can help with call-off orders where you buy and own the product, but we store it in our warehouse and 'call off' delivery to you, as and when you need it. Alternatively, you may place a standing order and have a predetermined quantity delivered at regular intervals. Contact us to discuss your options.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk

# All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|---|---|
|  50µl Tips for Qiagen® QIAgility<br>**Volume**: 50 µl<br>**Pack Size**: 960 Tips (10 Racks × 96 Tips)<br>Create PDF Data Sheet | 960 Tips (10 Racks × 96 Tips) | E1172-1418 |
|  50µl Tips for Qiagen® QIAgility (Case Size)<br>**Volume**: 50 µl<br>**Pack Size**: 4,800 Tips<br>(50 Racks × 96 Tips)<br>Create PDF Data Sheet | 4,800 Tips (50 Racks × 96 Tips) | E1172-1418-C |
|  200µl Tips for Qiagen® QIAgility<br>**Volume**: 200 µl<br>**Pack Size**: 960 Tips (10 Racks × 96 Tips)<br>Create PDF Data Sheet | 960 Tips (10 Racks × 96 Tips) | E1173-2418 |
|  200µl Tips for Qiagen® QIAgility (Case Size)<br>**Volume**: 200 µl<br>**Pack Size**: 4,800 Tips<br>(50 Racks × 96 Tips)<br>Create PDF Data Sheet | 4,800 Tips (50 Racks × 96 Tips) | E1173-2418-C |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk