# StarChill MCT, Purple/Pink



**Product Highlights**

- 24-well rack with lid, for use with: 0.5, 1.5 & 2.0 ml microcentrifuge / screw cap tubes
- Temperature colour change: 7 °C
- Estimated time at 4 °C: 3+ hours*
- Dimensions (WxDxH):
- 146 x 105 x 65 mm (without lid)
- 146 x 105 x 70 mm (with lid)
- StarChill Rack must be frozen for a minimum of 8 hours at -20°C before use
- Freeze the StarChill Rack upside down in order for the gel to freeze evenly

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

## General Data

| | |
|---|---|
| Art. No. | E2396-5634 |
| Dimensions (W × D × H) | 146 mm × 105 mm × 70 mm |
| Dimensions (W × D) | 146 mm × 105 mm |
| Autoclavable | No |
| Pack Size | 1 Piece (1 Box × 1 Piece) |
| Number of wells | 96 wells |
| Rack Type | Chiller Rack |
| Color | Pink (warm)/purple (cold) |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

# More informations about StarChill MCT, Purple/Pink

**StarChill high performance chiller racks**

Maintain a fridge environment on the benchtop! The centre section of StarChill racks changes colour when frozen. It then gradually turns back to its original colour as it warms up, thus giving an instant, visual-indication of the temperature.

*Time will vary, depending on ambient temperature. All StarChill racks must be frozen for a minimum of 8 hours at -20 °C before use. Not Autoclavable.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

# All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|---|---|
|  StarChill MCT, Purple/Pink<br>**Color**: Pink (warm)/purple (cold)<br>Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | E2396-5634 |
|  StarChill MCT, Green/Yellow<br>**Color**: Yellow (warm)/green (cold)<br>Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | E2396-5636 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de