200 µl XL Graduated Filter Tip / Wide Orifice, Racked (sterile)



Product Highlights

- Large orifice tip ideal for pipetting macromolecules and viscous samples
- Compatible with a wide range of pipettes
- Certified
- 1.0 mm inside orifice
- Supplied in convenient, press to open rack

Product Applications

Next generation sequencing. Cell culture. Plant sciences.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

General Data

Art. No.	E1011-7420
Length	91.1 mm
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume	200 µl
Pack Size	192 Tips (2 Racks × 96 Tips)
Tip Type	XL wide orifice
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

For use with brand + product (SL)	**Abimed:** BioPette-E, Discovery Comfort, LABMATE Optima **Biohit®:** eLine, ePet, mLine, Picus, Proline **BRAND®:** Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S **Capp®:** CappAero Comfort, CappAero Multi, CappAero Universal, Capp Bravo, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel **Eppendorf®:** Reference, Reference 2, Research, Research 3 neo, Research plus, Research pro, Xplorer, Xplorer Plus **Gilson®:** Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra **Hamilton®:** SoftGrip Pipette, SofTouch **HTL:** Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette **Nichiryo®:** 00-NP7, 50F, NLE, NPLO, NPX **Rainin®:** Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus **Sartorius®:** eLine, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta **Socorex®:** Acura 815 manual, Acura 825, Acura electro 925, Calibra digital **STARLAB®:** ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E **Thermo Fisher Scientific Inc.:** Finn Focus, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

More informations about 200 µl XL Graduated Filter Tip / Wide Orifice, Racked (sterile)

Protect delicate samples from shear force and reduce flow resistance with Starlab's wide orifice tips.

› Eliminate cell fragmentation caused by mechanical shearing

› Maintaining the integrity of high molecular weight DNA like gDNA or rDNA

› Pipetting of viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup,etc.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de