# 1000 µl XL Wide Orifice Tips



**Product Highlights**

- Large orifice tip ideal for pipetting macromolecules and viscous samples
- Compatible with a wide range of pipettes
- Certified
- 1.63 mm inside orifice

## Product Applications

Next generation sequencing. Cell culture. Plant sciences.
Eliminate cell fragmentation caused by mechanical shearing
Maintaining the integrity of high molecular weight DNA like gDNA or rDNA
Pipetting of viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup,etc.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

# General Data

| | |
|---|---|
| Art. No. | See variations |
| Filter tip | Non-filter |
| Volume | 1000 µl |
| Tip Type | XL wide orifice |
| Material | Polypropylene |
| DNA free | Yes |
| DNase free | Yes |
| PCR inhibitor free | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Graduations | Yes |
| Color | Natural |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

# More informations about 1000 µl XL Wide Orifice Tips

› Protect delicate samples from shear force and reduce flow resistance
› Eliminate cell fragmentation caused by mechanical shearing
› Maintaining the integrity of high molecular weight DNA like gDNA or rDNA
› Pipetting of viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup,etc.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

# All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|---|---|
|  1000 µl XL Graduated Tip / Wide Orifice, Loose (non-sterile)<br>**Sterile**: No<br>**Packaging Unit**: Bag<br>Create PDF Data Sheet | 1,000 Tips (1 Bag × 1000 Tips) | E1011-9000 |
|  1000 µl XL Graduated Tip / Wide Orifice, Racked (non-sterile)<br>**Sterile**: No<br>**Packaging Unit**: Rack<br>Create PDF Data Sheet | 960 Tips (10 Rack × 96 Tips) | E1011-9400 |
|  1000 µl XL Graduated Tip / Wide Orifice, Racked (sterile)<br>**Sterile**: Yes<br>**Packaging Unit**: Racks<br>Create PDF Data Sheet | 960 Tips (10 Racks × 96 Tips) | E1011-9410 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de