# ISO-Rack for 24 x 0.5 ml Microcentrifuge Tubes



**Product Highlights**

- ISO-Rack for 24 x 0.5 ml Tubes
- For use with Mixer HC Thermoblock S8012-0010

## Product Applications

Sample Preparation
Purification of plasmids, DNA, RNA
cDNA synthesis
Enzyme reactions (restriction digest, ligations, proteinase-K-digestion)
Denaturation of DNA, RNA, proteins
Labeling of DNA, RNA, proteins
Cultivation of bacteria in microcentrifuge tubes, centrifuge tubes and deepwell plates
Transformation of bacteria
Resuspension of pellets

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

# General Data

| | |
|---|---|
| Art. No. | S8012-0020 |
| Pack Size | 1 Piece (1 Box × 1 Piece) |
| Operating temperature | 4 to 35 °C |
| Sample capacity | 24 × 0.5 ml microcentrifuge tubes |
| Speed Range | 300 – 1,500 rpm |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

# More informations about ISO-Rack for 24 x 0.5 ml Microcentrifuge Tubes

ISO rack for use with 0.5ml tubes (Mixer HC Thermoblock S8012-0010). The Iso-Rack is ideal for time-critical applications and high throughput sample preparation. Prepare your samples in your tubes, using the Iso-Rack as your bench rack. You can then transfer all 24 tubes to the thermoblock at the same time, in seconds!

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de